



How to Use Stella with Pulse




Step 1

Assemble Stella with Pulse on a level surface using Stella cradle, close to a sink.




Step 2

Pre-clean instrument and devices thoroughly. Flush pre-clean solution through lumens.




Step 3

Push the 'ON' button on the Stella IQ and Stella Pulse, until you hear a sound. A Bluetooth connection will automatically occur. Follow the instructions on the Stella IQ.




Step 4

Add the pre-cleaned instrument and small parts tray into the inner instrument compartment of the Base.




Step 5

Connect the tube set to Pulse then connect the luer lock to the instrument channel. Place the tube mount on the tray cut out so the pick-up/filter is placed in the inner compartment.




Step 6

Make up 5 litres of Tristel Fuse for Instrument Disinfectant. Read Fuse for Instruments User Guide carefully.




Step 7

Pour the Tristel Fuse for Instruments disinfectant into the inner instrument compartment until the liquid overflows. Ensure you use all 5 litres of disinfectant. A beep confirms the start of the five minute cycle.




Step 8

Place the Lid onto the Base. This will force water to flow over the side walls of the instrument compartment and create a complete seal. The Stella IQ displays the disinfecting segment automatically counting down five minutes. Stella IQ directs Stella Pulse to pump disinfectant through the channel of the instrument for one minute.




Step 9

After 5 minutes Stella drains the disinfectant automatically. Stella Pulse pumps air through the instrument's channel to remove residual disinfectant. When the Stella Base is empty, there is an audible reminder that the cycle is complete.



Step 10

Press the OK button to confirm removal of the instrument and to acknowledge that the cycle is complete.



Step 11

Record the validation code in the Patient's Notes or the Audit Trail Record Book. Stella IQ will store all information ready for download to your computer. A new cycle can immediately be started by pressing the OK button again. To turn off, press the OFF button on Stella IQ.