

Stella Suite Bluetooth Installation User Guide



Version 1.2

© Copyright Reserved 2010 Stella Performance Ltd

Stella Suite



Bluetooth Installation User Guide



Table of Contents

1. Initiate Bluetooth on Stella5 IQ	3
2. Bluetooth Connections	4
3. Determine Port number.	6
4. Setup Stella Suite.....	7
5. Product Support Service.....	8

1. Initiate Bluetooth on Stella5 IQ

Connect Bluetooth

- 1.1. Operator press “ON” for duration to activate the Stella system

Note: This initial delay is to stop electrical interference and accidental contact to the “ON” button causing unintentional activation of the Stella system.

- 1.2. The display will light up and display the “STELLA” logo;

Graphic “Stella” Logo



- 1.3. The Logo will remain for 1 second.
During this time the operator must press the “ON” button a second time to activate the Bluetooth Download mode.
- 1.4. Sound - Beep.
- 1.5. The following graphics are displayed illustrating the Bluetooth download mode is ready.
Graphic Icon - Bluetooth mode ready



- 1.6. The Bluetooth Radio module is now active, ready to be linked to the PC.


Bluetooth link is connected to the PC




When the communication link is made with laptop or PC the Radio icon on the right hand side will be animated as displayed.

2. Bluetooth Connections

- 2.1. Select the “Start” button.
- 2.2. From the start menu select the “Control Panel” on the right of the start menu.



- 2.3. The control panel will open.
- 2.4. Double click on the Bluetooth icon.




2.5. Press the “Add” button under the Devices tab.

2.6. Tick the “My device is setup and ready to be found” tick box on the following screen.

2.7. Press the “Next” Button.

The computer will automatically start searching for all available Bluetooth devices.

If the Stella5 IQ is not found ensure the Stella5 IQ is still on and press the “Search Again” button.



2.8. Once the device is display either click on the “Stella” Icon and press the “Next” button or double click on the “Stella” Icon to install the Unit.



2.9. The computer will now request a pass key.

2.10. Select choose my own pass key.

2.11. Type in the box “1234”.

2.12. Press “Next” button.

2.13. The computer will install the Stella5 IQ connection.



2.14. Take note of the “Outgoing COM Port number”.

2.15. Press the “Finish” Button.

Outgoing COM Port is
COM4 in this case.




3. Determine Port number.

3.1. Follow steps 2.1 to 2.4 in this document.

3.2. Click on the Stella5 IQ Icon and then click properties.

3.3. Select the Services or Devices tab and the COM Port is displayed in the List.



4. Setup Stella Suite


Bluetooth

This allows the user to manually configure the Stella Suite Bluetooth communication port if required.

Tick on the “Use Bluetooth” option for wireless radio link download.

Consult your Bluetooth instruction manual on how to create a Bluetooth communication link from the PC to the Stella5 IQ.

A COM port will be issued when searching for the Stella5 IQ Bluetooth link. Please select from the list of COM ports the Bluetooth COM port number.



- 4.1. Open Stella Suite and select “File – Configure”
- 4.2. Select the “Bluetooth” tab.
- 4.3. Select the COM port the Stella IQ is setup for.

In the example above we have used “COM4”

- 4.4. Tick the “Use Bluetooth” box.
- 4.5. Click “OK”.
- 4.6. To download press “the Download” button.



5. Product Support Service

Please note if you have questions regarding the operation of this product please read this installation and user guide carefully before making further enquiries.

For more information on this product or other products please contact;

Stella Performance Limited Offices

New Zealand

PO Box 4276
Mt Maunganui 3149,
New Zealand.
T: +64 (7) 575 3730
F: +64 (7) 575 3760
E: info@stella-performance.com
W: www.stella-performance.com

China

Shanghai Stella Medical Equipment Co.Ltd
Rm, 19, 24th Floor, 227 Zun Yi Road
Changning District
Shanghai
China.
T: +86 (21) 325 02858
F: +86 (21) 325 02858
E: info@stella-performance.com
W: www.stella-performance.com

United Kingdom

Unit 1B, Lynx Business Park,
Fordham Road,
Snailwell, Cambridgeshire,
CB8 7NY, UK.
T: +44 (0)1638 721500
F: +44 (0)1638 721911
E: info@stella-performance.com
W: www.stella-performance.com

Manual Versions

Date	Version	Author	Changes	Check	Approve
05/05/09	1.2	B Kent	Revise Contact details, Front cover page, Meta Normal-Roman Font, Removed reference to Vista, amended graphics	JB	AH
04/02/09	1.1	DROGERS	Original Version	B Kent	AW

STADC207